

Technical Data Sheet

YAO-YUAN ECG Gel — 140 g Tube

Manufacturer: Guangzhou Yao-Yuan Industrial Co., Ltd.
Brand: YAO-YUAN
Address: Building 4, No. 3 Huashan Road, Shilou Town, Panyu District, Guangzhou 511447, Guangzhou, Guangdong, China

Description:

YAO-YUAN ECG Gel is a conductive, water-based gel designed for use in electrocardiography (ECG) procedures. It ensures reliable signal transmission between the skin and electrodes, while being gentle to the skin, odorless, and easy to clean.

Key Features:

Odorless	Transparent / Colorless
Water-soluble; easy to remove	Non-abrasive; leaves no residue
Supplied in 140 g tube	Compatible with common ECG electrodes

Technical Details:

Formulation	Water-based, non-ionic
pH	Skin-friendly
Viscosity	Optimized for controlled spread and contact
Salt content	Salt-free, non-corrosive to electrodes
Storage	Keep tightly closed, in a cool, dry place; avoid freezing and sunlight
Shelf life	As indicated on packaging (batch specific)

Packaging & Labeling:

Primary container: 140 g tube. Label includes lot/batch number, manufacturing date, expiry date, country of origin, and brand name.

Contact:

Guangzhou Yao-Yuan Industrial Co., Ltd.
Building 4, No. 3 Huashan Road, Shilou Town, Panyu District,
Guangzhou 511447, Guangzhou, Guangdong, China
Email: info@yaoyuangel.com / sales@yaoyuangel.com